



PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

PRODUCT CATALOG

Tissue Expanders,
Breast Implants and Sizers

EFFECTIVE JULY 2023

WARNING:

- Breast implants are not considered lifetime devices. The longer people have them, the greater the chances are that they will develop complications, some of which will require more surgery.
- Breast implants have been associated with the development of a cancer of the immune system called breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. Some patients have died from BIA-ALCL.
- Patients receiving breast implants have reported a variety of systemic symptoms such as joint pain, muscle aches, confusion, chronic fatigue, autoimmune diseases and others. Individual patient risk for developing these symptoms has not been well established. Some patients report complete resolution of symptoms when the implants are removed without replacement.

The sale and distribution of Mentor Breast Implant Devices are restricted to users and/or user facilities that provide information to patients about the risks and benefits of the device prior to its use in the form and manner specified in approved labeling to be provided by Mentor Worldwide LLC.

CONTENTS

MemoryGel® Breast Implants	4-12
and Resterilizable Smooth Gel Sizers	
Smooth Round Moderate Plus Profile BOOST Breast Implant	
Smooth Round Moderate High Profile BOOST Breast Implant	
Smooth Round High Profile BOOST Breast Implant	
Smooth Round Moderate Plus Profile Xtra Breast Implant	
Smooth Round Moderate High Profile Xtra Breast Implant	
Smooth Round High Profile Xtra Breast Implant	
Smooth Round Moderate Classic Profile Breast Implant	
Smooth Round Moderate Plus Profile Breast Implant	
Smooth Round High Profile Breast Implant	
Smooth Round Ultra High Profile Breast Implant	
SILTEX® Microtextured Round Mod Plus Xtra Breast Implant	
SILTEX® Microtextured Round High Profile Xtra Breast Implant	
SILTEX® Microtextured Round Moderate Classic Profile Breast Implant	
SILTEX® Microtextured Round Moderate Plus Profile Breast Implant	
SILTEX® Microtextured Round High Profile Breast Implant	
SILTEX® Microtextured Round Ultra High Profile Breast Implant	
MemoryShape® SILTEX® Microtextured Breast Implants.....	13-15
and Resterilizable Gel Sizers	
MM Profile, Medium Height/Moderate Projection Breast Implant	
MM+ Profile, Medium Height/Moderate Plus Projection Breast Implant	
MH Profile, Medium Height/High Projection Breast Implant	
TM+ Profile, Tall Height/Moderate Plus Projection Breast Implant	
LM+ Profile, Low Height/Moderate Plus Projection Breast Implant	
TH Profile, Tall Height/High Projection Breast Implant	
Saline Breast Implants.....	16-17
and Sterile Saline Breast Sizers	
Smooth Round Moderate Profile Breast Implant	
Smooth Round Moderate Plus Profile Breast Implant	
Smooth Round High Profile Breast Implant	
Round Moderate Profile Breast Sizer	
Round Moderate Plus Profile Breast Sizer	
Round High Profile Breast Sizer	
SPECTRUM® Adjustable Saline Breast Implants	18
Smooth Round SPECTRUM® Saline Breast Implant	
Tissue Expanders.....	19-22
ARTOURA® Breast Tissue Expanders	
Smooth and SILTEX® High Profile Breast Tissue Expanders	
Smooth and SILTEX® Ultra High Profile Breast Tissue Expanders	
CPX® 4 Breast Tissue Expanders	
Smooth and SILTEX® Low Height Breast Tissue Expanders	
Smooth and SILTEX® Medium Height Breast Tissue Expanders	
Smooth and SILTEX® Tall Height Breast Tissue Expanders	
Smooth Tissue Expanders	
Smooth Round Expander	
Smooth Elliptical Expander	
Smooth Crescent Expander	
Smooth Rectangle Expander	
Accessories.....	23
MENTOR® Volume Sizing System	
Aseptic Transfer System	
Injection Port Detector	
Injection Domes	
Miscellaneous	
MENTOR® Promise Protection Plan	24
and MENTOR® Deflation Protection Program	
Important Safety Information	25

MemoryGel® BOOST™ Smooth Breast Implants



Smooth Round Moderate Plus Profile BOOST

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
130 cc	8.8 cm	3.1 cm	SMPB130	RSZMPB130S
160 cc	9.3 cm	3.3 cm	SMPB160	RSZMPB160S
185 cc	9.8 cm	3.3 cm	SMPB185	RSZMPB185S
210 cc	10.3 cm	3.5 cm	SMPB210	RSZMPB210S
235 cc	10.7 cm	3.6 cm	SMPB235	RSZMPB235S
265 cc	11.2 cm	3.6 cm	SMPB265	RSZMPB265S
290 cc	11.4 cm	3.9 cm	SMPB290	RSZMPB290S
320 cc	11.6 cm	4.1 cm	SMPB320	RSZMPB320S
340 cc	12.0 cm	4.1 cm	SMPB340	RSZMPB340S
390 cc	12.7 cm	4.1 cm	SMPB390	RSZMPB390S
430 cc	12.9 cm	4.4 cm	SMPB430	RSZMPB430S
450 cc	13.2 cm	4.5 cm	SMPB450	RSZMPB450S
480 cc	13.4 cm	4.6 cm	SMPB480	RSZMPB480S
530 cc	13.9 cm	4.7 cm	SMPB530	RSZMPB530S
545 cc	14.1 cm	4.8 cm	SMPB545	RSZMPB545S
590 cc	14.3 cm	5.0 cm	SMPB590	RSZMPB590S
630 cc	14.8 cm	5.0 cm	SMPB630	RSZMPB630S
675 cc	15.1 cm	5.0 cm	SMPB675	RSZMPB675S
730 cc	15.6 cm	5.1 cm	SMPB730	RSZMPB730S

MemoryGel® BOOST™ Smooth Breast Implants



Smooth Round High Profile BOOST

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
200 cc	9.0 cm	4.4 cm	SHPB200	RSZHPB200S
235 cc	9.5 cm	4.7 cm	SHPB235	RSZHPB235S
260 cc	9.8 cm	4.8 cm	SHPB260	RSZHPB260S
285 cc	10.3 cm	4.9 cm	SHPB285	RSZHPB285S
325 cc	10.5 cm	5.3 cm	SHPB325	RSZHPB325S
370 cc	11.0 cm	5.4 cm	SHPB370	RSZHPB370S
405 cc	11.2 cm	5.6 cm	SHPB405	RSZHPB405S
440 cc	11.6 cm	5.8 cm	SHPB440	RSZHPB440S
480 cc	12.2 cm	5.8 cm	SHPB480	RSZHPB480S
520 cc	12.3 cm	6.0 cm	SHPB520	RSZHPB520S
545 cc	12.5 cm	6.1 cm	SHPB545	RSZHPB545S
585 cc	12.8 cm	6.3 cm	SHPB585	RSZHPB585S
635 cc	13.2 cm	6.4 cm	SHPB635	RSZHPB635S
685 cc	13.7 cm	6.4 cm	SHPB685	RSZHPB685S
740 cc	14.1 cm	6.6 cm	SHPB740	RSZHPB740S
775 cc	14.4 cm	6.6 cm	SHPB775	RSZHPB775S

Smooth Round Moderate High Profile BOOST



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
150 cc	8.8 cm	3.7 cm	SMHB150	RSZMHB150S
170 cc	9.2 cm	3.9 cm	SMHB170	RSZMHB170S
190 cc	9.5 cm	3.9 cm	SMHB190	RSZMHB190S
210 cc	9.8 cm	4.2 cm	SMHB210	RSZMHB210S
230 cc	10.1 cm	4.3 cm	SMHB230	RSZMHB230S
260 cc	10.5 cm	4.3 cm	SMHB260	RSZMHB255S
280 cc	10.7 cm	4.4 cm	SMHB280	RSZMHB280S
310 cc	11.1 cm	4.6 cm	SMHB310	RSZMHB310S
330 cc	11.3 cm	4.7 cm	SMHB330	RSZMHB330S
355 cc	11.6 cm	4.8 cm	SMHB355	RSZMHB350S
375 cc	11.8 cm	4.9 cm	SMHB375	RSZMHB375S
400 cc	12.1 cm	5.0 cm	SMHB400	RSZMHB400S
435 cc	12.4 cm	5.2 cm	SMHB435	RSZMHB430S
470 cc	12.8 cm	5.2 cm	SMHB470	RSZMHB470S
500 cc	13.2 cm	5.3 cm	SMHB500	RSZMHB500S
545 cc	13.5 cm	5.5 cm	SMHB545	RSZMHB545S
580 cc	13.8 cm	5.5 cm	SMHB580	RSZMHB580S
635 cc	14.2 cm	5.6 cm	SMHB635	RSZMHB630S
700 cc	14.8 cm	5.8 cm	SMHB700	RSZMHB700S
775 cc	15.3 cm	6.0 cm	SMHB775	RSZMHB775S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

*Shipped sterile. Resterilizable up to 10 times.

MemoryGel® Xtra Smooth Breast Implants and Resterilizable Gel Sizers

Smooth Round Moderate Plus Profile Xtra

Volume	Diameter	Projection	Catalog #
115 cc	8.3 cm	3.1 cm	SMPX115
130 cc	8.9 cm	3.1 cm	SMPX130
160 cc	9.5 cm	3.3 cm	SMPX160
190 cc	10.0 cm	3.4 cm	SMPX190
215 cc	10.4 cm	3.6 cm	SMPX215
240 cc	10.8 cm	3.7 cm	SMPX240
270 cc	11.4 cm	3.8 cm	SMPX270
295 cc	11.5 cm	4.0 cm	SMPX295
325 cc	11.9 cm	4.1 cm	SMPX325
350 cc	12.3 cm	4.1 cm	SMPX350
370 cc	12.6 cm	4.1 cm	SMPX370
405 cc	12.7 cm	4.4 cm	SMPX405
440 cc	13.1 cm	4.5 cm	SMPX440
465 cc	13.3 cm	4.7 cm	SMPX465
490 cc	13.6 cm	4.7 cm	SMPX490
525 cc	13.9 cm	4.7 cm	SMPX525
545 cc	14.0 cm	4.9 cm	SMPX545
560 cc	14.3 cm	4.8 cm	SMPX560
605 cc	14.5 cm	5.1 cm	SMPX605
630 cc	14.7 cm	5.1 cm	SMPX630
645 cc	14.9 cm	5.1 cm	SMPX645
685 cc	15.3 cm	5.2 cm	SMPX685
755 cc	15.7 cm	5.4 cm	SMPX755



Resterilizable Gel Sizer*

RSZMPX115S
RSZMPX130S
RSZMPX160S
RSZMPX190S
RSZMPX215S
RSZMPX240S
RSZMPX270S
RSZMPX295S
RSZMPX325S
RSZMPX350S
RSZMPX370S
RSZMPX405S
RSZMPX440S
RSZMPX465S
RSZMPX490S
RSZMPX525S
RSZMPX545S
RSZMPX560S
RSZMPX605S
RSZMPX630S
RSZMPX645S
RSZMPX685S
RSZMPX755S

Smooth Round Moderate High Profile Xtra

Volume	Diameter	Projection	Catalog #
130 cc	8.4 cm	3.5 cm	SMHX130
150 cc	8.8 cm	3.7 cm	SMHX150
170 cc	9.1 cm	3.8 cm	SMHX170
185 cc	9.5 cm	3.9 cm	SMHX185
210 cc	9.8 cm	4.1 cm	SMHX210
230 cc	10.1 cm	4.2 cm	SMHX230
255 cc	10.5 cm	4.3 cm	SMHX255
280 cc	10.8 cm	4.4 cm	SMHX280
310 cc	11.1 cm	4.6 cm	SMHX310
330 cc	11.4 cm	4.7 cm	SMHX330
350 cc	11.6 cm	4.8 cm	SMHX350
375 cc	11.9 cm	4.9 cm	SMHX375
400 cc	12.1 cm	4.9 cm	SMHX400
430 cc	12.5 cm	5.1 cm	SMHX430
470 cc	12.8 cm	5.2 cm	SMHX470
500 cc	13.1 cm	5.3 cm	SMHX500
545 cc	13.4 cm	5.5 cm	SMHX545
580 cc	13.8 cm	5.5 cm	SMHX580
630 cc	14.2 cm	5.7 cm	SMHX630
700 cc	14.6 cm	5.9 cm	SMHX700
775 cc	15.2 cm	6.0 cm	SMHX775



Resterilizable Gel Sizer*

RSZSMHX130S
RSZSMHX150S
RSZSMHX170S
RSZSMHX185S
RSZSMHX210S
RSZSMHX230S
RSZSMHX255S
RSZSMHX280S
RSZSMHX310S
RSZSMHX330S
RSZSMHX350S
RSZSMHX375S
RSZSMHX400S
RSZSMHX430S
RSZSMHX470S
RSZSMHX500S
RSZSMHX545S
RSZSMHX580S
RSZSMHX630S
RSZSMHX700S
RSZSMHX775S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

{6}

*Shipped sterile. Resterilizable up to 10 times.

MemoryGel® Xtra Smooth Breast Implants and Resterilizable Gel Sizers

Smooth Round High Profile Xtra

Volume	Diameter	Projection	Catalog #
150 cc	8.3 cm	4.1 cm	SHPX150
180 cc	8.8 cm	4.4 cm	SHPX180
200 cc	9.2 cm	4.5 cm	SHPX200
235 cc	9.5 cm	4.9 cm	SHPX235
260 cc	9.9 cm	5.0 cm	SHPX260
285 cc	10.3 cm	5.0 cm	SHPX285
335 cc	10.7 cm	5.4 cm	SHPX335
355 cc	11.0 cm	5.5 cm	SHPX355
380 cc	11.3 cm	5.5 cm	SHPX380
415 cc	11.5 cm	5.8 cm	SHPX415
450 cc	11.9 cm	5.9 cm	SHPX450
465 cc	12.1 cm	5.8 cm	SHPX465
490 cc	12.5 cm	5.8 cm	SHPX490
535 cc	12.7 cm	6.0 cm	SHPX535
560 cc	12.7 cm	6.3 cm	SHPX560
595 cc	13.1 cm	6.3 cm	SHPX595
650 cc	13.5 cm	6.5 cm	SHPX650
700 cc	14.1 cm	6.5 cm	SHPX700
755 cc	14.4 cm	6.7 cm	SHPX755
790 cc	14.8 cm	6.7 cm	SHPX790



Resterilizable Gel Sizer*

RSZSHPX150S
RSZSHPX180S
RSZSHPX200S
RSZSHPX235S
RSZSHPX260S
RSZSHPX285S
RSZSHPX335S
RSZSHPX355S
RSZSHPX380S
RSZSHPX415S
RSZSHPX450S
RSZSHPX465S
RSZSHPX490S
RSZSHPX535S
RSZSHPX560S
RSZSHPX595S
RSZSHPX650S
RSZSHPX700S
RSZSHPX755S
RSZSHPX790S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

{7}

*Shipped sterile. Resterilizable up to 10 times.

MemoryGel® Smooth Breast Implants and Resterilizable Gel Sizers



Smooth Round Moderate Classic Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
130 cc	9.6 cm	2.4 cm	3507130 MC	RSZ7130MCS
150 cc	10.3 cm	2.5 cm	3507150 MC	RSZ7150MCS
170 cc	10.5 cm	2.6 cm	3507170 MC	RSZ7170MCS
190 cc	11.0 cm	2.7 cm	3507190 MC	RSZ7190MCS
215 cc	11.6 cm	2.8 cm	3507215 MC	RSZ7215MCS
235 cc	11.8 cm	2.9 cm	3507235 MC	RSZ7235MCS
255 cc	12.2 cm	2.9 cm	3507255 MC	RSZ7255MCS
275 cc	12.5 cm	3.0 cm	3507275 MC	RSZ7275MCS
300 cc	12.8 cm	3.1 cm	3507300 MC	RSZ7300MCS
320 cc	13.0 cm	3.2 cm	3507320 MC	RSZ7320MCS
340 cc	13.4 cm	3.2 cm	3507340 MC	RSZ7340MCS
360 cc	13.6 cm	3.3 cm	3507360 MC	RSZ7360MCS
385 cc	13.9 cm	3.3 cm	3507385 MC	RSZ7385MCS
405 cc	14.2 cm	3.4 cm	3507405 MC	RSZ7405MCS
425 cc	14.6 cm	3.4 cm	3507425 MC	RSZ7425MCS
445 cc	14.8 cm	3.5 cm	3507445 MC	RSZ7445MCS
470 cc	15.0 cm	3.5 cm	3507470 MC	RSZ7470MCS
490 cc	15.2 cm	3.7 cm	3507490 MC	RSZ7490MCS
510 cc	15.4 cm	3.7 cm	3507510 MC	RSZ7510MCS
555 cc	15.7 cm	3.8 cm	3507555 MC	RSZ7555MCS
595 cc	16.1 cm	3.8 cm	3507595 MC	RSZ7595MCS
640 cc	16.5 cm	3.9 cm	3507640 MC	RSZ7640MCS
680 cc	17.0 cm	3.9 cm	3507680 MC	RSZ7680MCS
750 cc	17.5 cm	4.0 cm	3507750 MC	RSZ7750MCS
800 cc	18.0 cm	4.1 cm	3507800 MC	RSZ7800MCS



Smooth Round Moderate Plus Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8.2 cm	2.7 cm	3501001 BC	RSZ1001S
125 cc	8.9 cm	2.8 cm	3501251 BC	RSZ1251S
150 cc	9.5 cm	2.9 cm	3501501 BC	RSZ1501S
175 cc	10.0 cm	3.1 cm	3501751 BC	RSZ1751S
200 cc	10.5 cm	3.2 cm	3502001 BC	RSZ2001S
225 cc	10.9 cm	3.3 cm	3502251 BC	RSZ2251S
250 cc	11.3 cm	3.4 cm	3502501 BC	RSZ2501S
275 cc	11.7 cm	3.5 cm	3502751 BC	RSZ2751S
300 cc	12.0 cm	3.6 cm	3503001 BC	RSZ3001S
325 cc	12.3 cm	3.8 cm	3503251 BC	RSZ3251S
350 cc	12.5 cm	3.9 cm	3503501 BC	RSZ3501S
375 cc	12.8 cm	4.0 cm	3503751 BC	RSZ3751S
400 cc	13.1 cm	4.0 cm	3504001 BC	RSZ4001S
425 cc	13.5 cm	4.1 cm	3504251 BC	RSZ4251S
450 cc	13.6 cm	4.2 cm	3504501 BC	RSZ4501S
475 cc	14.0 cm	4.2 cm	3504751 BC	RSZ4751S
500 cc	14.1 cm	4.3 cm	3505001 BC	RSZ5001S
525 cc	14.5 cm	4.4 cm	3505251 BC	RSZ5251S
550 cc	14.6 cm	4.5 cm	3505501 BC	RSZ5501S
575 cc	14.9 cm	4.6 cm	3505751 BC	RSZ5751S
600 cc	15.0 cm	4.6 cm	3506001 BC	RSZ6001S
650 cc	15.4 cm	4.8 cm	3506501 BC	RSZ6501S
700 cc	15.8 cm	4.9 cm	3507001 BC	RSZ7001S
750 cc	16.3 cm	4.9 cm	3507501 BC	RSZ7501S
800 cc	16.5 cm	5.1 cm	3508001 BC	RSZ8001S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

*Shipped sterile. Resterilizable up to 10 times.

{8}

MemoryGel® Smooth Breast Implants and Resterilizable Gel Sizers



Smooth Round High Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8.3 cm	3.5 cm	3501254 BC	RSZ1254S
150 cc	8.8 cm	3.7 cm	3501504 BC	RSZ1504S
175 cc	9.3 cm	3.9 cm	3501754 BC	RSZ1754S
200 cc	9.7 cm	4.0 cm	3502004 BC	RSZ2004S
225 cc	10.1 cm	4.2 cm	3502254 BC	RSZ2254S
250 cc	10.5 cm	4.3 cm	3502504 BC	RSZ2504S
275 cc	10.8 cm	4.4 cm	3502754 BC	RSZ2754S
300 cc	11.1 cm	4.5 cm	3503004 BC	RSZ3004S
325 cc	11.4 cm	4.6 cm	3503254 BC	RSZ3254S
350 cc	11.7 cm	4.8 cm	3503504 BC	RSZ3504S
375 cc	12.0 cm	4.8 cm	3503754 BC	RSZ3754S
400 cc	12.2 cm	5.0 cm	3504004 BC	RSZ4004S
425 cc	12.5 cm	5.0 cm	3504254 BC	RSZ4254S
450 cc	12.8 cm	5.1 cm	3504504 BC	RSZ4504S
475 cc	12.9 cm	5.3 cm	3504754 BC	RSZ4754S
500 cc	13.2 cm	5.3 cm	3505004 BC	RSZ5004S
550 cc	13.6 cm	5.5 cm	3505504 BC	RSZ5504S
600 cc	14.0 cm	5.6 cm	3506004 BC	RSZ6004S
650 cc	14.4 cm	5.7 cm	3506504 BC	RSZ6504S
700 cc	14.8 cm	5.8 cm	3507004 BC	RSZ7004S
750 cc	15.1 cm	5.9 cm	3507504 BC	RSZ7504S
800 cc	15.5 cm	6.0 cm	3508004 BC	RSZ8004S



Smooth Round Ultra High Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
135 cc	7.8 cm	4.1 cm	3505135 BC	RSZ5135S
160 cc	8.2 cm	4.3 cm	3505160 BC	RSZ5160S
185 cc	8.4 cm	4.4 cm	3505185 BC	RSZ5185S
215 cc	8.7 cm	4.5 cm	3505215 BC	RSZ5215S
240 cc	9.0 cm	4.7 cm	3505240 BC	RSZ5240S
270 cc	9.3 cm	4.8 cm	3505270 BC	RSZ5270S
295 cc	9.7 cm	5.0 cm	3505295 BC	RSZ5295S
320 cc	9.8 cm	5.1 cm	3505320 BC	RSZ5320S
350 cc	10.1 cm	5.2 cm	3505350 BC	RSZ5350S
375 cc	10.4 cm	5.3 cm	3505375 BC	RSZ5375S
400 cc	10.6 cm	5.4 cm	3505400 BC	RSZ5400S
430 cc	10.9 cm	5.6 cm	3505430 BC	RSZ5430S
455 cc	11.2 cm	5.7 cm	3505455 BC	RSZ5455S
480 cc	11.4 cm	5.8 cm	3505480 BC	RSZ5480S
535 cc	12.0 cm	6.1 cm	3505535 BC	RSZ5535S
590 cc	12.5 cm	6.3 cm	3505590 BC	RSZ5590S
650 cc	12.9 cm	6.4 cm	3505650 BC	RSZ5650S
700 cc	13.3 cm	6.5 cm	3505700 BC	RSZ5700S
750 cc	13.7 cm	6.6 cm	3505750 BC	RSZ5750S
800 cc	14.1 cm	6.7 cm	3505800 BC	RSZ5800S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

*Shipped sterile. Resterilizable up to 10 times.

{9}

MemoryGel® SILTEX® Microtextured Xtra Breast Implants and Resterilizable Gel Sizers

SILTEX® Microtextured Round Moderate Plus Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
115 cc	8.4 cm	3.2 cm	TMPX115	RSZMPX115S
130 cc	9.1 cm	3.1 cm	TMPX130	RSZMPX130S
160 cc	9.5 cm	3.3 cm	TMPX160	RSZMPX160S
190 cc	10.0 cm	3.5 cm	TMPX190	RSZMPX190S
215 cc	10.5 cm	3.6 cm	TMPX215	RSZMPX215S
240 cc	10.9 cm	3.8 cm	TMPX240	RSZMPX240S
270 cc	11.4 cm	3.9 cm	TMPX270	RSZMPX270S
295 cc	11.6 cm	4.1 cm	TMPX295	RSZMPX295S
325 cc	11.9 cm	4.2 cm	TMPX325	RSZMPX325S
350 cc	12.2 cm	4.3 cm	TMPX350	RSZMPX350S
370 cc	12.7 cm	4.3 cm	TMPX370	RSZMPX370S
405 cc	12.8 cm	4.8 cm	TMPX405	RSZMPX405S
440 cc	13.1 cm	4.8 cm	TMPX440	RSZMPX440S
465 cc	13.4 cm	4.8 cm	TMPX465	RSZMPX465S
490 cc	13.6 cm	4.8 cm	TMPX490	RSZMPX490S
525 cc	13.7 cm	5.1 cm	TMPX525	RSZMPX525S
545 cc	14.1 cm	5.0 cm	TMPX545	RSZMPX545S
560 cc	14.3 cm	5.0 cm	TMPX560	RSZMPX560S
605 cc	14.4 cm	5.3 cm	TMPX605	RSZMPX605S
630 cc	14.6 cm	5.4 cm	TMPX630	RSZMPX630S
645 cc	15.0 cm	5.3 cm	TMPX645	RSZMPX645S
685 cc	15.4 cm	5.3 cm	TMPX685	RSZMPX685S
755 cc	15.7 cm	5.6 cm	TMPX755	RSZMPX755S

SILTEX® Microtextured Round High Profile Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer†
150 cc	8.4 cm	4.2 cm	THPX150	RSZTHPX150S
175 cc	9.0 cm	4.4 cm	THPX175	RSZTHPX175S
200 cc	9.4 cm	4.5 cm	THPX200	RSZTHPX200S
230 cc	9.8 cm	4.8 cm	THPX230	RSZTHPX230S
255 cc	10.1 cm	4.9 cm	THPX255	RSZTHPX255S
285 cc	10.5 cm	5.1 cm	THPX285	RSZTHPX285S
325 cc	10.8 cm	5.4 cm	THPX325	RSZTHPX325S
340 cc	11.1 cm	5.4 cm	THPX340	RSZTHPX340S
365 cc	11.4 cm	5.5 cm	THPX365	RSZTHPX365S
405 cc	11.6 cm	5.9 cm	THPX405	RSZTHPX405S
425 cc	11.9 cm	5.8 cm	THPX425	RSZTHPX425S
455 cc	12.1 cm	6.0 cm	THPX455	RSZTHPX455S
470 cc	12.3 cm	5.9 cm	THPX470	RSZTHPX470S
515 cc	12.6 cm	6.2 cm	THPX515	RSZTHPX515S
545 cc	12.8 cm	6.3 cm	THPX545	RSZTHPX545S
570 cc	13.2 cm	6.3 cm	THPX570	RSZTHPX570S
620 cc	13.4 cm	6.5 cm	THPX620	RSZTHPX620S
680 cc	14.0 cm	6.6 cm	THPX680	RSZTHPX680S
725 cc	14.5 cm	6.7 cm	THPX725	RSZTHPX725S
765 cc	14.7 cm	6.7 cm	THPX765	RSZTHPX765S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

*Shipped sterile. Resterilizable up to 10 times.

MemoryGel® SILTEX® Microtextured Breast Implants and Resterilizable Gel Sizers

SILTEX® Microtextured Round Moderate Classic Profile



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
130 cc	9.6 cm	2.5 cm	3541307MC	RSZ7130MCS
150 cc	10.2 cm	2.6 cm	3541507MC	RSZ7150MCS
170 cc	10.6 cm	2.7 cm	3541707MC	RSZ7170MCS
190 cc	11.0 cm	2.8 cm	3541907MC	RSZ7190MCS
215 cc	11.5 cm	2.9 cm	3542157MC	RSZ7215MCS
235 cc	11.7 cm	3.0 cm	3542357MC	RSZ7235MCS
255 cc	12.1 cm	3.1 cm	3542557MC	RSZ7255MCS
275 cc	12.3 cm	3.2 cm	3542757MC	RSZ7275MCS
300 cc	12.6 cm	3.3 cm	3543007MC	RSZ7300MCS
320 cc	12.9 cm	3.3 cm	3543207MC	RSZ7320MCS
340 cc	13.2 cm	3.4 cm	3543407MC	RSZ7340MCS
360 cc	13.4 cm	3.4 cm	3543607MC	RSZ7360MCS
385 cc	13.7 cm	3.5 cm	3543857MC	RSZ7385MCS
405 cc	14.0 cm	3.6 cm	3544057MC	RSZ7405MCS
425 cc	14.3 cm	3.6 cm	3544257MC	RSZ7425MCS
445 cc	14.5 cm	3.7 cm	3544457MC	RSZ7445MCS
470 cc	14.7 cm	3.7 cm	3544707MC	RSZ7470MCS
490 cc	14.9 cm	3.8 cm	3544907MC	RSZ7490MCS
510 cc	15.1 cm	3.8 cm	3545107MC	RSZ7510MCS
555 cc	15.4 cm	4.0 cm	3545557MC	RSZ7555MCS
595 cc	15.8 cm	4.1 cm	3545957MC	RSZ7595MCS
640 cc	16.3 cm	4.2 cm	3546407MC	RSZ7640MCS
680 cc	16.7 cm	4.2 cm	3546807MC	RSZ7680MCS
750 cc	17.2 cm	4.3 cm	3547507MC	RSZ7750MCS
800 cc	17.6 cm	4.4 cm	3548007MC	RSZ7800MCS

SILTEX® Microtextured Round Moderate Plus Profile



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8.1 cm	2.7 cm	3541001	RSZ1001S
125 cc	8.8 cm	2.9 cm	3541251	RSZ1251S
150 cc	9.4 cm	3.0 cm	3541501	RSZ1501S
175 cc	10.0 cm	3.2 cm	3541751	RSZ1751S
200 cc	10.5 cm	3.3 cm	3542001	RSZ2001S
225 cc	10.9 cm	3.5 cm	3542251	RSZ2251S
250 cc	11.3 cm	3.6 cm	3542501	RSZ2501S
275 cc	11.7 cm	3.7 cm	3542751	RSZ2751S
300 cc	12.0 cm	3.7 cm	3543001	RSZ3001S
325 cc	12.3 cm	3.8 cm	3543251	RSZ3251S
350 cc	12.6 cm	3.8 cm	3543501	RSZ3501S
375 cc	12.9 cm	3.9 cm	3543751	RSZ3751S
400 cc	13.2 cm	4.0 cm	3544001	RSZ4001S
425 cc	13.5 cm	4.2 cm	3544251	RSZ4251S
450 cc	13.7 cm	4.1 cm	3544501	RSZ4501S
475 cc	13.9 cm	4.4 cm	3544751	RSZ4751S
500 cc	14.1 cm	4.2 cm	3545001	RSZ5001S
525 cc	14.4 cm	4.6 cm	3545251	RSZ5251S
550 cc	14.4 cm	4.4 cm	3545501	RSZ5501S
575 cc	14.9 cm	4.7 cm	3545751	RSZ5751S
600 cc	14.7 cm	4.5 cm	3546001	RSZ6001S
650 cc	15.5 cm	4.9 cm	3546501	RSZ6501S
700 cc	15.7 cm	4.8 cm	3547001	RSZ7001S
750 cc	16.2 cm	5.1 cm	3547501	RSZ7501S
800 cc	16.6 cm	5.0 cm	3548001	RSZ8001S

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

*Shipped sterile. Resterilizable up to 10 times.

MemoryGel® SILTEX® Microtextured Breast Implants
and Resterilizable Gel Sizers



SILTEX® Microtextured Round High Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8.4 cm	3.6 cm	3544125	RSZ1254S
150 cc	8.9 cm	3.8 cm	3544150	RSZ1504S
175 cc	9.4 cm	4.0 cm	3544175	RSZ1754S
200 cc	9.9 cm	4.1 cm	3544200	RSZ2004S
225 cc	10.2 cm	4.3 cm	3544225	RSZ2254S
250 cc	10.5 cm	4.5 cm	3544250	RSZ2504S
275 cc	10.9 cm	4.6 cm	3544275	RSZ2754S
300 cc	11.1 cm	4.7 cm	3544300	RSZ3004S
325 cc	11.5 cm	4.8 cm	3544325	RSZ3254S
350 cc	11.7 cm	4.9 cm	3544350	RSZ3504S
375 cc	12.0 cm	5.0 cm	3544375	RSZ3754S
400 cc	12.3 cm	5.1 cm	3544400	RSZ4004S
425 cc	12.5 cm	5.2 cm	3544425	RSZ4254S
450 cc	12.7 cm	5.2 cm	3544450	RSZ4504S
475 cc	12.9 cm	5.4 cm	3544475	RSZ4754S
500 cc	13.2 cm	5.4 cm	3544500	RSZ5004S
550 cc	13.5 cm	5.6 cm	3544550	RSZ5504S
600 cc	14.0 cm	5.7 cm	3544600	RSZ6004S
650 cc	14.3 cm	5.8 cm	3544650	RSZ6504S
700 cc	14.7 cm	6.0 cm	3544700	RSZ7004S
750 cc	15.1 cm	6.0 cm	3544750	RSZ7504S
800 cc	15.4 cm	6.3 cm	3544800	RSZ8004S

SILTEX® Microtextured Round Ultra High Profile

Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
135 cc	8.0 cm	4.3 cm	3545135	RSZ5135S
160 cc	8.4 cm	4.5 cm	3545160	RSZ5160S
185 cc	8.6 cm	4.6 cm	3545185	RSZ5185S
215 cc	8.9 cm	4.7 cm	3545215	RSZ5215S
240 cc	9.2 cm	4.9 cm	3545240	RSZ5240S
270 cc	9.5 cm	5.0 cm	3545270	RSZ5270S
295 cc	9.9 cm	5.2 cm	3545295	RSZ5295S
320 cc	10.0 cm	5.3 cm	3545320	RSZ5320S
350 cc	10.3 cm	5.4 cm	3545350	RSZ5350S
375 cc	10.5 cm	5.5 cm	3545375	RSZ5375S
400 cc	10.8 cm	5.6 cm	3545400	RSZ5400S
430 cc	11.1 cm	5.8 cm	3545430	RSZ5430S
455 cc	11.3 cm	5.9 cm	3545455	RSZ5455S
480 cc	11.6 cm	6.0 cm	3545480	RSZ5480S
535 cc	12.1 cm	6.3 cm	3545535	RSZ5535S
590 cc	12.6 cm	6.5 cm	3545590	RSZ5590S
650 cc	13.0 cm	6.6 cm	3545650	RSZ5650S
700 cc	13.4 cm	6.7 cm	3545700	RSZ5700S

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

*Shipped sterile.
Resterilizable up to 10 times.

MemoryShape® SILTEX® Microtextured Breast Implants
and Resterilizable Gel Sizers



MM Profile, Medium Height/Moderate Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
120 cc	9.0 cm	8.5 cm	3.3 cm	3540908	RSZ0908S
135 cc	9.5 cm	8.9 cm	3.5 cm	3540958	RSZ0958S
155 cc	10.0 cm	9.4 cm	3.7 cm	3541008	RSZ1008S
180 cc	10.5 cm	9.9 cm	3.8 cm	3541058	RSZ1058S
215 cc	11.0 cm	10.3 cm	3.9 cm	3541108	RSZ1108S
245 cc	11.5 cm	10.8 cm	4.0 cm	3541158	RSZ1158S
280 cc	12.0 cm	11.3 cm	4.2 cm	3541208	RSZ1208S
315 cc	12.5 cm	11.8 cm	4.4 cm	3541258	RSZ1258S
355 cc	13.0 cm	12.2 cm	4.6 cm	3541308	RSZ1308S
395 cc	13.5 cm	12.7 cm	4.7 cm	3541358	RSZ1358S
440 cc	14.0 cm	13.2 cm	4.9 cm	3541408	RSZ1408S
480 cc	14.5 cm	13.6 cm	5.0 cm	3541458	RSZ1458S
530 cc	15.0 cm	14.1 cm	5.2 cm	3541508	RSZ1508S
640 cc	16.0 cm	15.0 cm	5.6 cm	3541608	RSZ1608S
775 cc	17.0 cm	16.0 cm	5.9 cm	3541708	RSZ1708S

MM+ Profile, Medium Height/
Moderate Plus Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
140 cc	9.0 cm	8.5 cm	3.8 cm	3340905	RSZ0905S
165 cc	9.5 cm	8.9 cm	4.0 cm	3340955	RSZ0955S
195 cc	10.0 cm	9.4 cm	4.2 cm	334-005	RSZ1005S
225 cc	10.5 cm	9.9 cm	4.4 cm	3341055	RSZ1055S
255 cc	11.0 cm	10.3 cm	4.7 cm	3341105	RSZ1105S
295 cc	11.5 cm	10.8 cm	4.9 cm	3341155	RSZ1155S
330 cc	12.0 cm	11.3 cm	5.1 cm	3341205	RSZ1205S
375 cc	12.5 cm	11.8 cm	5.3 cm	3341255	RSZ1255S
420 cc	13.0 cm	12.2 cm	5.5 cm	3341305	RSZ1305S
475 cc	13.5 cm	12.7 cm	5.7 cm	3341355	RSZ1355S
525 cc	14.0 cm	13.2 cm	5.9 cm	3341405	RSZ1405S
585 cc	14.5 cm	13.6 cm	6.1 cm	3341455	RSZ1455S
650 cc	15.0 cm	14.1 cm	6.3 cm	3341505	RSZ1505S

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

*Shipped sterile.
Resterilizable up to 10 times.

MemoryShape® SILTEX® Microtextured Breast Implants
and Resterilizable Gel Sizers



MH Profile, Medium Height/High Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
165 cc	9.0 cm	8.5 cm	4.6 cm	3340902	RSZ0902S
195 cc	9.5 cm	8.9 cm	4.8 cm	3340952	RSZ0952S
225 cc	10.0 cm	9.4 cm	5.1 cm	3341002	RSZ1002S
260 cc	10.5 cm	9.9 cm	5.3 cm	3341052	RSZ1052S
300 cc	11.0 cm	10.3 cm	5.6 cm	3341102	RSZ1102S
345 cc	11.5 cm	10.8 cm	5.8 cm	3341152	RSZ1152S
390 cc	12.0 cm	11.3 cm	6.0 cm	3341202	RSZ1202S
440 cc	12.5 cm	11.8 cm	6.2 cm	3341252	RSZ1252S
495 cc	13.0 cm	12.2 cm	6.5 cm	3341302	RSZ1302S
555 cc	13.5 cm	12.7 cm	6.7 cm	3341352	RSZ1352S
620 cc	14.0 cm	13.2 cm	6.9 cm	3341402	RSZ1402S
685 cc	14.5 cm	13.6 cm	7.1 cm	3341452	RSZ1452S

MemoryShape® SILTEX® Microtextured Breast Implants
and Resterilizable Gel Sizers



TM+ Profile, Tall Height/Moderate Plus Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
145 cc	9.0 cm	9.4 cm	3.8 cm	3340909	RSZ0909S
175 cc	9.5 cm	9.9 cm	4.0 cm	3340959	RSZ0959S
205 cc	10.0 cm	10.4 cm	4.2 cm	3341009	RSZ1009S
235 cc	10.5 cm	10.9 cm	4.4 cm	3341059	RSZ1059S
270 cc	11.0 cm	11.5 cm	4.7 cm	3341109	RSZ1109S
305 cc	11.5 cm	12.0 cm	4.9 cm	3341159	RSZ1159S
350 cc	12.0 cm	12.5 cm	5.1 cm	3341209	RSZ1209S
395 cc	12.5 cm	13.0 cm	5.3 cm	3341259	RSZ1259S
445 cc	13.0 cm	13.5 cm	5.5 cm	3341309	RSZ1309S
495 cc	13.5 cm	14.1 cm	5.7 cm	3341359	RSZ1359S
555 cc	14.0 cm	14.6 cm	5.9 cm	334409	RSZ1409S
615 cc	14.5 cm	15.1 cm	6.1 cm	3341459	RSZ1459S
680 cc	15.0 cm	15.6 cm	6.3 cm	3341509	RSZ1509S

TH Profile, Tall Height/High Projection



Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
180 cc	9.0 cm	9.4 cm	4.6 cm	3340904	RSZ4090S
215 cc	9.5 cm	9.9 cm	4.8 cm	3340954	RSZ4095S
250 cc	10.0 cm	10.4 cm	5.1 cm	3341004	RSZ4100S
290 cc	10.5 cm	10.9 cm	5.3 cm	3341054	RSZ4105S
330 cc	11.0 cm	11.5 cm	5.6 cm	3341104	RSZ4110S
380 cc	11.5 cm	12.0 cm	5.8 cm	3341154	RSZ4115S
430 cc	12.0 cm	12.5 cm	6.0 cm	3341204	RSZ4120S
485 cc	12.5 cm	13.0 cm	6.2 cm	3341254	RSZ4125S
545 cc	13.0 cm	13.5 cm	6.5 cm	3341304	RSZ4130S
610 cc	13.5 cm	14.1 cm	6.7 cm	3341354	RSZ4135S
680 cc	14.0 cm	14.6 cm	6.9 cm	3341404	RSZ4140S
755 cc	14.5 cm	15.1 cm	7.1 cm	3341454	RSZ4145S

LM+ Profile, Low Height/
Moderate Plus Projection



Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
125 cc	9.0 cm	8.0 cm	3.8 cm	3340907	RSZ0907S
145 cc	9.5 cm	8.4 cm	4.0 cm	3340957	RSZ0957S
170 cc	10.0 cm	8.8 cm	4.2 cm	3341007	RSZ1007S
195 cc	10.5 cm	9.3 cm	4.4 cm	3341057	RSZ1057S
225 cc	11.0 cm	9.7 cm	4.7 cm	3341107	RSZ1107S
255 cc	11.5 cm	10.2 cm	4.9 cm	3341157	RSZ1157S
290 cc	12.0 cm	10.6 cm	5.1 cm	3341207	RSZ1207S
330 cc	12.5 cm	11.1 cm	5.3 cm	3341257	RSZ1257S
370 cc	13.0 cm	11.5 cm	5.5 cm	3341307	RSZ1307S
415 cc	13.5 cm	12.0 cm	5.7 cm	3341357	RSZ1357S
465 cc	14.0 cm	12.4 cm	5.9 cm	3341407	RSZ1407S
515 cc	14.5 cm	12.8 cm	6.1 cm	3341457	RSZ1457S
570 cc	15.0 cm	13.3 cm	6.3 cm	3341507	RSZ1507S
690 cc	16.0 cm	14.2 cm	6.8 cm	3341607	RSZ1607S

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

*Shipped sterile.
Resterilizable up to 10 times.

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

*Shipped sterile.
Resterilizable up to 10 times.

SALINE BREAST IMPLANTS

Smooth Round Moderate Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer†
Volume	Diameter	Proj.	Vol.	Diam.	Proj.		
125 cc	9.5 cm	3.0 cm	150 cc	9.1 cm	3.4 cm	3501610	351125SZ
150 cc	10.0 cm	3.1 cm	175 cc	9.6 cm	3.7 cm	3501615	
175 cc	10.6 cm	3.3 cm	200 cc	10.0 cm	3.9 cm	3501620	351175SZ
200 cc	11.0 cm	3.4 cm	225 cc	10.4 cm	3.9 cm	3501625	
225 cc	11.5 cm	3.5 cm	250 cc	10.9 cm	4.0 cm	3501630	351225SZ
250 cc	11.9 cm	3.6 cm	275 cc	11.2 cm	4.1 cm	3501635	
275 cc	12.3 cm	3.7 cm	300 cc	11.7 cm	4.2 cm	3501640	351275SZ
300 cc	12.6 cm	3.7 cm	325 cc	12.1 cm	4.3 cm	3501645	
325 cc	13.0 cm	3.8 cm	375 cc	12.3 cm	4.5 cm	3501650	351325SZ
350 cc	13.3 cm	3.9 cm	400 cc	12.6 cm	4.7 cm	3501655	
375 cc	13.6 cm	4.0 cm	425 cc	12.8 cm	4.8 cm	3501660	351375SZ
425 cc	14.2 cm	4.1 cm	475 cc	13.5 cm	4.9 cm	3501670	351425SZ
475 cc	14.8 cm	4.2 cm	525 cc	14.2 cm	4.9 cm	3501680	351475SZ
525 cc	15.0 cm	4.2 cm	575 cc	14.5 cm	5.1 cm	3501685	
575 cc	15.0 cm	4.5 cm	625 cc	14.8 cm	5.1 cm	3501690	
625 cc	15.2 cm	4.6 cm	700 cc	14.9 cm	5.6 cm	3501695	
700 cc	15.6 cm	4.9 cm	775 cc	14.9 cm	5.8 cm	3501697	

†Single Use Only

Smooth Round Moderate Plus Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer†
Volume	Diameter	Proj.	Vol.	Diam.	Proj.		
175 cc	9.5 cm	3.5 cm	210 cc	9.3 cm	4.2 cm	3502175	3512175SZ
200 cc	10.0 cm	3.7 cm	240 cc	9.7 cm	4.4 cm	3502200	3512200SZ
225 cc	10.4 cm	3.8 cm	270 cc	10.1 cm	4.6 cm	3502225	3512225SZ
250 cc	10.8 cm	4.0 cm	300 cc	10.5 cm	4.7 cm	3502250	3512250SZ
275 cc	11.0 cm	4.1 cm	330 cc	10.8 cm	4.9 cm	3502275	3512275SZ
300 cc	11.5 cm	4.3 cm	360 cc	11.2 cm	5.0 cm	3502300	3512300SZ
325 cc	11.9 cm	4.4 cm	390 cc	11.5 cm	5.2 cm	3502325	3512325SZ
350 cc	12.1 cm	4.5 cm	420 cc	11.7 cm	5.3 cm	3502350	3512350SZ
375 cc	12.3 cm	4.6 cm	450 cc	12.0 cm	5.4 cm	3502375	3512375SZ
400 cc	12.6 cm	4.7 cm	480 cc	12.3 cm	5.5 cm	3502400	3512400SZ
425 cc	12.9 cm	4.8 cm	510 cc	12.5 cm	5.6 cm	3502425	3512425SZ
450 cc	13.0 cm	4.9 cm	540 cc	12.8 cm	5.7 cm	3502450	3512450SZ
475 cc	13.3 cm	5.0 cm	570 cc	13.0 cm	5.8 cm	3502475	3512475SZ
500 cc	13.6 cm	5.1 cm	600 cc	13.2 cm	5.9 cm	3502500	3512500SZ
550 cc	14.0 cm	5.3 cm	660 cc	13.7 cm	6.1 cm	3502550	3512550SZ
600 cc	14.5 cm	5.5 cm	720 cc	14.1 cm	6.3 cm	3502600	3512600SZ
650 cc	14.8 cm	5.6 cm	780 cc	14.4 cm	6.5 cm	3502650	3512650SZ
700 cc	15.2 cm	5.8 cm	840 cc	14.8 cm	6.6 cm	3502700	3512700SZ
750 cc	15.6 cm	5.9 cm	900 cc	15.1 cm	6.8 cm	3502750	3512750SZ
800 cc	16.1 cm	6.0 cm	960 cc	15.5 cm	7.0 cm	3502800	3512800SZ

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

†Single Use Only

SALINE BREAST IMPLANTS

Smooth Round High Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer†
Volume	Diameter	Proj.	Vol.	Diam.	Proj.		
170 cc	9.0 cm	4.0 cm	200 cc	8.8 cm	4.9 cm	3503170	
190 cc	9.3 cm	4.1 cm	225 cc	9.1 cm	5.1 cm	3503190	3513190SZ
210 cc	9.6 cm	4.2 cm	250 cc	9.5 cm	5.3 cm	3503210	
230 cc	10.0 cm	4.3 cm	275 cc	9.8 cm	5.5 cm	3503230	3513230SZ
250 cc	10.2 cm	4.5 cm	300 cc	10.0 cm	5.6 cm	3503250	
270 cc	10.4 cm	4.6 cm	325 cc	10.2 cm	5.8 cm	3503270	3513270SZ
290 cc	10.8 cm	4.7 cm	350 cc	10.5 cm	5.9 cm	3503290	3513290SZ
310 cc	11.0 cm	4.8 cm	375 cc	10.7 cm	6.0 cm	3503310	3513310SZ
330 cc	11.3 cm	4.8 cm	400 cc	11.0 cm	6.2 cm	3503330	3513330SZ
380 cc	11.7 cm	5.2 cm	450 cc	11.4 cm	6.4 cm	3503380	3513380SZ
420 cc	12.0 cm	5.4 cm	500 cc	11.7 cm	6.7 cm	3503420	3513420SZ
460 cc	12.4 cm	5.5 cm	550 cc	12.1 cm	6.9 cm	3503460	
500 cc	12.8 cm	5.6 cm	600 cc	12.4 cm	7.1 cm	3503500	3513500SZ
560 cc	13.4 cm	5.7 cm	675 cc	13.1 cm	7.3 cm	3503560	
630 cc	13.8 cm	5.9 cm	750 cc	13.4 cm	7.4 cm	3503630	3513630SZ

†Single Use Only

Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.

SPECTRUM® POST-OPERATIVELY ADJUSTABLE SALINE BREAST IMPLANTS

Smooth Round SPECTRUM®


NOMINAL FILL			MAXIMUM FILL	Catalog #
Volume	Diameter	Projection	Volume	
125 cc	9.4 cm	2.7 cm	150 cc	3501410
175 cc	10.3 cm	3.1 cm	210 cc	3501420
225 cc	10.8 cm	3.5 cm	270 cc	3501430
275 cc	11.6 cm	3.8 cm	330 cc	3501440
325 cc	12.4 cm	4.0 cm	390 cc	3501450
375 cc	13.5 cm	4.1 cm	450 cc	3501460
425 cc	14.0 cm	4.2 cm	510 cc	3501470
475 cc	14.3 cm	4.4 cm	570 cc	3501480
525 cc	15.0 cm	4.4 cm	630 cc	3501485
575 cc	14.7 cm	4.8 cm	690 cc	3501490




BREAST TISSUE EXPANDERS

ARTOURA® PLUS Smooth and SILTEX® Breast Tissue Expanders High Profile and Ultra High Profile

All expanders include a CENTERSCOPE® Magnetic Detection Device and 21-Gauge Winged Infusion Set



High Profile				PLUS Smooth Surface 6 Suture Tabs	SILTEX® Texture Surface 3 Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
225 cc	10.0 cm	10.0 cm	5.2 cm	SDC100H	TEXP100RH
300 cc	11.0 cm	11.0 cm	6.0 cm	SDC110H	TEXP110RH
375 cc	12.0 cm	12.0 cm	6.4 cm	SDC120H	TEXP120RH
475 cc	13.0 cm	13.0 cm	6.8 cm	SDC130H	TEXP130RH
500 cc	13.5 cm	13.5 cm	6.6 cm	SDC135H	TEXP135RH
600 cc	14.0 cm	14.0 cm	7.1 cm	SDC140H	TEXP140RH
750 cc	15.0 cm	15.0 cm	7.7 cm	SDC150H	TEXP150RH
850 cc	15.5 cm	15.5 cm	7.8 cm	SDC155H	TEXP155RH



Ultra High Profile				PLUS Smooth Surface 6 Suture Tabs	SILTEX® Texture Surface 3 Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
350 cc	10.0 cm	10.0 cm	7.0 cm	SDC100UH	TEXP100RUH
455 cc	11.0 cm	11.0 cm	7.6 cm	SDC110UH	TEXP110RUH
535 cc	12.0 cm	12.0 cm	8.0 cm	SDC120UH	TEXP120RUH
650 cc	13.0 cm	13.0 cm	8.2 cm	SDC130UH	TEXP130RUH
700 cc	13.5 cm	13.5 cm	8.3 cm	SDC135UH	TEXP135RUH
850 cc	14.0 cm	14.0 cm	8.9 cm	SDC140UH	TEXP140RUH

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.


PLEASE NOTE: CPX® and ARTOURA™ Breast Tissue Expanders are not MRI compatible.

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.


BREAST TISSUE EXPANDERS

CPX®4 PLUS Smooth Breast Tissue Expander with Suture Tabs


All expanders include a CENTERSCOPE® Magnetic Detection Device and 21-Gauge Winged Infusion Set



Low Height				PLUS Smooth Surface 6 Suture Tabs
Volume	Width (cm)	Height (cm)	Projection	Catalog #
250 cc	11.4 cm	8.1 cm	6.1 cm	SCPX114LH
350 cc	12.7 cm	9.4 cm	6.5 cm	SCPX127LH
450 cc	14.0 cm	10.2 cm	7.1 cm	SCPX140LH
550 cc	15.0 cm	10.9 cm	7.4 cm	SCPX150LH
650 cc	15.7 cm	11.2 cm	7.9 cm	SCPX157LH
750 cc	16.5 cm	11.9 cm	8.1 cm	SCPX165LH



Medium Height				PLUS Smooth Surface 6 Suture Tabs
Volume	Width (cm)	Height (cm)	Projection	Catalog #
275 cc	10.7 cm	9.3 cm	6.2 cm	SCPX107MH
350 cc	11.7 cm	10.0 cm	6.6 cm	SCPX117MH
450 cc	12.7 cm	10.8 cm	7.0 cm	SCPX127MH
550 cc	13.5 cm	11.7 cm	7.4 cm	SCPX135MH
650 cc	14.6 cm	12.6 cm	7.6 cm	SCPX146MH
800 cc	15.6 cm	13.3 cm	8.0 cm	SCPX156MH



Tall Height				PLUS Smooth Surface 6 Suture Tabs
Volume	Width (cm)	Height (cm)	Projection	Catalog #
250 cc	10.1 cm	10.7 cm	5.6 cm	SCPX101TH
350 cc	11.3 cm	11.8 cm	6.0 cm	SCPX113TH
450 cc	12.3 cm	12.9 cm	6.5 cm	SCPX123TH
550 cc	13.2 cm	13.8 cm	6.9 cm	SCPX132TH
650 cc	14.0 cm	14.6 cm	7.3 cm	SCPX140TH
750 cc	14.6 cm	15.3 cm	7.6 cm	SCPX146TH
850 cc	15.4 cm	15.9 cm	7.9 cm	SCPX154TH

PLEASE NOTE: CPX® Breast Tissue Expanders are not MRI compatible.

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

BREAST TISSUE EXPANDERS

CPX®4 SILTEX® Breast Tissue Expander & CPX®4 SILTEX® Breast Tissue Expander with Suture Tabs


All expanders include a CENTERSCOPE® Magnetic Detection Device and 21-Gauge Winged Infusion Set



Low Height				SILTEX® Texture Surface 3 Suture Tabs	SILTEX® Texture Surface No Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
250 cc	11.4 cm	8.1 cm	6.1 cm	3549111	3548111
350 cc	12.7 cm	9.4 cm	6.5 cm	3549112	3548112
450 cc	14.0 cm	10.2 cm	7.1 cm	3549113	3548113
550 cc	15.0 cm	10.9 cm	7.4 cm	3549114	3548114
650 cc	15.7 cm	11.2 cm	7.9 cm	3549115	3548115
750 cc	16.5 cm	11.9 cm	8.1 cm	3549116	3548116



Medium Height				SILTEX® Texture Surface 3 Suture Tabs	SILTEX® Texture Surface No Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
275 cc	10.7 cm	9.3 cm	6.2 cm	3549211	3548211
350 cc	11.7 cm	10.0 cm	6.6 cm	3549212	3548212
450 cc	12.7 cm	10.8 cm	7.0 cm	3549213	3548213
550 cc	13.5 cm	11.7 cm	7.4 cm	3549214	3548214
650 cc	14.6 cm	12.6 cm	7.6 cm	3549215	3548215
800 cc	15.6 cm	13.3 cm	8.0 cm	3549216	3548216



Tall Height				SILTEX® Texture Surface 3 Suture Tabs	SILTEX® Texture Surface No Suture Tabs
Volume	Width	Height	Projection	Catalog #	Catalog #
250 cc	10.1 cm	10.7 cm	5.6 cm	3549311	3548311
350 cc	11.3 cm	11.8 cm	6.0 cm	3549312	3548312
450 cc	12.3 cm	12.9 cm	6.5 cm	3549313	3548313
550 cc	13.2 cm	13.8 cm	6.9 cm	3549314	3548314
650 cc	14.0 cm	14.6 cm	7.3 cm	3549315	3548315
750 cc	14.6 cm	15.3 cm	7.6 cm	3549316	3548316
850 cc	15.4 cm	15.9 cm	7.9 cm	3549317	3548317
1445 cc	14.9 cm	16.0 cm	11.6 cm	3549336	

PLEASE NOTE: CPX® and ARTOURA™ Breast Tissue Expanders are not MRI compatible.

Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

TISSUE EXPANDERS

Smooth Tissue Expanders

Smooth Round

Device Volume	Diameter	Projection	Catalog #
400 cc	11.3 cm	5.9 cm	3504305M
550 cc	13.4 cm	5.8 cm	3504307M
700 cc	14.6 cm	6.7 cm	3504304M
850 cc	15.0 cm	7.2 cm	3504308M



Smooth Elliptical

Device Volume	Length	Width	Nominal Projection	TEMPORARY OVEREXPANSION		Catalog #
				Proj.	Max. Vol.	
25 cc	6.0 cm	4.0 cm	2.3 cm	4.5 cm	75 cc	3505301M
50 cc	8.0 cm	5.0 cm	2.5 cm	5.3 cm	150 cc	3505302M
75 cc	10.0 cm	6.0 cm	2.1 cm	6.0 cm	225 cc	3505303M
150 cc	12.0 cm	7.0 cm	3.2 cm	7.2 cm	450 cc	3505304M
250 cc	14.0 cm	8.0 cm	4.5 cm	9.0 cm	750 cc	3505305M
125 cc	15.0 cm	5.0 cm	3.3 cm	6.5 cm	375 cc	3505306M




Smooth Crescent

Device Volume	Length	Width	Projection	Catalog #
50 cc	9.0 cm	5.7 cm	2.7 cm	3504312M
100 cc	11.0 cm	6.9 cm	3.7 cm	3504313M
200 cc	14.6 cm	7.5 cm	4.3 cm	3504314M
300 cc	16.5 cm	10.3 cm	4.7 cm	3504315M
400 cc	18.2 cm	11.4 cm	5.2 cm	3504316M
500 cc	19.6 cm	12.4 cm	5.5 cm	3504317M



Smooth Rectangle

Device Volume	Length	Width	Projection	Catalog #
50 cc	6.9 cm	3.2 cm	3.2 cm	3504309M
100 cc	6.6 cm	5.1 cm	5.6 cm	3504301M
250 cc	10.4 cm	7.3 cm	5.6 cm	3504302M
400 cc	10.6 cm	9.3 cm	6.6 cm	3504303M
550 cc	13.7 cm	7.1 cm	7.7 cm	3504310M
700 cc	15.0 cm	8.0 cm	7.9 cm	3504311M



Individual implant dimensions may vary slightly in products of this type.
Not all units will conform exactly to the dimensions noted above.

ACCESSORIES

MENTOR® Volume Sizing System


Carton contains four sterile kits. Each kit contains the following:

Description	Catalog #
MENTOR® Volume Sizing System - Display Case + Volume Sizers	ESZKIT
Replacement Units	Catalog #
Display Case	XXX
50 cc Partial External "Gap Filling" Sizer*	ESZ50
75 cc Partial External "Gap Filling" Sizer*	ESZ75
125 cc Partial External "Gap Filling" Sizer*	ESZ125
150 cc External Implant Sizer*	ESZ150
250 cc External Implant Sizer*	ESZ250
350 cc External Implant Sizer*	ESZ350
450 cc External Implant Sizer*	ESZ450
550 cc External Implant Sizer*	ESZ550
650 cc External Implant Sizer*	ESZ650

Aseptic Transfer System

Carton contains six sterile kits. Each kit contains the following:

Description	Catalog #
Aseptic Transfer System	3508600
Includes:	
• 50 cc syringe (sterile, disposable)	
• Piercing spike	
• 48" tubing (saline container to syringe)	
• Dual check valve	
• 10" tubing (syringe to implant fill tube)	
• Luer connector (attaches to implant fill tube)	



Injection Port Detector

Description	Catalog #
CENTERSCOPE® Magnetic Detection Device	3504402

Injection Domes

Description	Catalog #
Dome Pack Set	350DOMPK

Miscellaneous

Description	Catalog #
Standard Fill Tube 8"	600975001
Long Fill Tube 24"	600975002
Check Valve	200464001
21-Gauge Winged Infusion Set	7B3050

*A customer is able to purchase replacement units only following the purchase of the MENTOR® Volume Sizing System (ESZKIT).

PLEASE NOTE: Individual implant dimensions may vary slightly in products of this type.
Not all units conform exactly to the dimensions noted above.

The MENTOR® | Promise Protection Plan The Best and Most Comprehensive Implant Warranty¹

MENTOR® Breast Implants are backed by the industry's most comprehensive warranty. MENTOR® Promise—with FREE and automatic comprehensive coverage—brings the peace of mind that our breast implants are protected for years to come.

The MENTOR® Deflation Protection Program² Mentor Family of Tissue Expanders

MENTOR® is dedicated to advancing plastic surgery by listening to surgeons and customers and responding to the ever changing needs of healthcare providers in the US. The MENTOR® Deflation Protection Program is backed by our Replacement Policy for CPX® and ARTOURA™ family of Tissue Expanders in the case of deflation of a MENTOR® Breast Tissue Expander in a patient.

For more information about the MENTOR® Promise Protection Plan and the MENTOR® Deflation Protection Program please visit www.breastimplantsbymentor.com/breast-implants/warranties

References:

1. Allergan Warranty http://www.natrellesurgeon.com/Content/assets/Natrelle_Warranty%20Brochure_DFU. Accessed 04/12/2018. Sientra Product replacement policy and warranty . http://www.sientra.com/Content/pdfs/Warranty_Terms_And_Conditions.pdf. Accessed 04/12/2018. Sientra CapCon Warranty. http://www.sientra.com/Content/pdfs/CapCon_Care_Program_Terms_and_Conditions.pdf. Accessed 04/12/2018. Sientra mailing to US Plastic Surgeons
2. Mentor Replacement Policy for the CPX® and Artoura™ Family of Tissue Expanders

For further information, talk to your Mentor sales representative or visit mentorwwllc.com

IMPORTANT SAFETY INFORMATION:

MENTOR® MemoryGel® Breast Implants, MENTOR® MemoryShape® Breast Implants, and MENTOR® Saline-filled Breast Implants are indicated for breast augmentation in women (at least 22 years old for MemoryGel® Implants and MemoryShape® Implants, and 18 years old for Saline Implants) or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body, with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions, or who are currently pregnant or nursing.

Breast implants are not lifetime devices and breast implantation may not be a one-time surgery. The chance of developing complications increases over time. The most common complications for breast augmentation and reconstruction with MemoryGel® Implants include any reoperation, capsular contracture, and implant removal with or without replacement. The most common complications with MemoryShape® Implants for breast augmentation include reoperation for any reason, implant removal with or without replacement, and ptosis. The most common complications with MemoryShape® Implants for breast reconstruction include reoperation for any reason, implant removal with or without replacement, and capsular contracture. A lower risk of complication is rupture. The health consequences of a ruptured silicone gel breast implant have not been fully established. MRI screenings are recommended three years after initial implant surgery and then every two years after to detect silent rupture. Breast implants are also associated with the risk of breast implant-associated anaplastic large cell lymphoma (BIA-ALCL), an uncommon type of lymphoma. An individual's risk of developing BIA-ALCL with MENTOR® Breast Implants is low based on the incidence of worldwide cases. The most common complications with MENTOR® Saline-filled Implants include reoperation, implant removal, capsular contracture, breast pain, and implant deflation.

For MemoryGel® Implants, patients should receive a copy of Important Information for Augmentation Patients about MENTOR® MemoryGel® Breast Implants or Important Information for Reconstruction Patients about MENTOR® MemoryGel® Breast Implants. For MemoryShape® Implants, patients should receive a copy of Patient Educational Brochure - Breast Augmentation with MENTOR® MemoryShape® Breast Implants or Patient Educational Brochure - Breast Reconstruction with MENTOR® MemoryShape® Breast Implants, and a copy of Quick Facts about Breast Augmentation & Reconstruction with MENTOR® MemoryShape® Breast Implants. For MENTOR® Saline-filled Implants, patients should receive a copy of Saline-Filled Breast Implants: Making an Informed Decision. Your patient needs to read and understand the information regarding the risks and benefits of breast implants, with an opportunity to consult with you prior to deciding on surgery.

The ARTOURA® Breast Tissue Expander or CONTOUR PROFILE® Breast Tissue Expander can be utilized for breast reconstruction after mastectomy, correction of an underdeveloped breast, scar revision, and tissue defect procedures. The expander is intended for temporary subcutaneous or submuscular implantation and is not intended for use beyond six months. ARTOURA® Breast Tissue Expander and CONTOUR PROFILE® Tissue Expanders include magnetic injection domes, which contain a rare earth permanent magnet, and are NOT MRI compatible. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. Do not use the ARTOURA® Tissue Expander nor CONTOUR PROFILE® Tissue Expander in patients where an MRI may be needed. DO NOT use the ARTOURA® Breast Tissue Expander and CONTOUR PROFILE® Tissue Expander in patients that have a previously implanted device such as pacemakers, drug infusion devices, artificial sensing devices, etc. that could be affected by a magnetic field. Mentor has not tested the in vivo effects of radiation therapy with ARTOURA® Breast Tissue Expander and CONTOUR PROFILE® Expander devices and cannot warrant the safety of such use. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas, where severe surgical reduction of the area has previously been performed and where steroids are used in the surgical pocket.

For detailed indications, contraindications, warnings, and precautions associated with the use of all MENTOR® Implantable Devices, which include MENTOR® Saline-filled Implants, MemoryGel® Implants, MemoryShape® Implants, ARTOURA® Expanders, and CONTOUR PROFILE® Expanders, please refer to the Instructions for Use (IFU) provided with each product or visit www.mentorwwllc.com.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, warnings and precautions.

Caution: US law restricts this device to sale by or on the order of a physician.



PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Mentor Worldwide, LLC

Irvine, CA 92618, USA

Customer Service: (800) 235-5731 | Fax: (866) 225-2873

© Mentor Worldwide LLC 2023 020827-230208